

HEADER BAGS

Medical Grade Header Bags Ensure Aseptic Presentation for Medical Device, Surgical Kits and Pharmaceutical Applications Requiring ETO or Radiation Sterilization.



Techni-Vent Header Bags are made of DuPont™ Medical Grade Tyvek® for strength and bacterial barrier, ensuring clean visual transfer of coating which confirms pouch integrity and peelability with no fiber tear or particulate. They are made of puncture and abrasion resistant film and laminates under strict ISO 13485 production standards. Our header bags can be manufactured to accommodate irregular and/or high-profile items.

www.technipaq.com

975 Lutter Drive
Crystal Lake, IL 60014

Phone: 815.477.1800
Fax: 815.477.0777
info@technipaq.com

•DUPONT•
Tyvek
Authorized converter of
medical grade Tyvek®



LAMINATING

COATING

PRINTING

SLITTING

CONVERTING



MEDICAL
DEVICE



DIAGNOSTIC



LIFE
SCIENCE



PHARMACEUTICAL

TO MEET YOUR SPECIFIC REQUIREMENTS

We offer a variety of films, laminates, and foil structures in a wide range of sizes.

CAPABILITIES

- Stock sizes available for immediate shipment – with Tyvek® vents
- Small or large quantity runs
- Unprinted or printed pouches in up to eight colors
- Custom pouches available up to 48"x99"

MATERIALS

- Nylon, polyethylene, or high puncture resistance options
- Specialty films
- DuPont™ Tyvek®
- Peelable laminates
- Foil laminates
- Clear high barrier

STERILIZATION OPTIONS

- Radiation
- Autoclave
- ETO

*We do
what others
don't.*

*We customize
where others
won't.*

STOCKED SIZES

Featuring 1073B DuPont™ Tyvek® with 4mil LLDPE

8" x 12" w/3.75"

10" x 15" w/3.75"

12" x 16" w/3.75"

13" x 15" w/3.75"

16" x 22" w/3.75"

19" x 30" w/3.75"



**DESIGN YOUR OWN
POUCH IN MINUTES**
technipaq.com/design

«DUPONT»
Tyvek
Authorized converter of
medical grade Tyvek®



www.technipaq.com

975 Lutter Drive
Crystal Lake, IL 60014

Phone: 815.477.1800
Fax: 815.477.0777
info@technipaq.com